

SuProne Plus™ System Attachments

Attachments for the durable positioning device.

The SuProne Plus System provides a number of different treatment options with three unique attachments. Each device attaches to the VersaBoard allowing for 360° treatment angle.

RT-7664 Carbon Fiber 45° Tilt Attachment

- Attaches to SuProne Plus Head & Neck Board
- May be positioned 5 degrees to 45 degrees
- Easy to locate and read degree location
- Stable base for maximum reproducibility
- Locks down to most couches
- Attaches to VersaBoard™

RT-7029 Carbon Fiber Prone Attachment (also available in acrylic - RT-5029)

- Attaches to SuProne Plus Head & Neck Board
- May be positioned -20 degrees or +20 degrees
- Includes comfort cushion for patient
- Easy to locate and read degree location
- Stable base for maximum reproducibility
- Locks down to most couches
- Attaches to VersaBoard™

RT-5049 Lateral Attachment

- Attaches to SuProne Plus™ Head and Neck Board
- Three standard positioning heights included (2cm, 5cm and 7cm)
- Can be used prone
- Locks down to most couches
- Attaches to VersaBoard™

RT- 7045 The VersaBoard™ Adapter

Allows you to attach the 45° and Prone to the VersaBoard™.


*To order, or for more information call **Bionix Radiation Therapy** at 1.800.624.6649 or visit us on the web at www.BionixRT.com.*



**45° attachment with
SuProne Plus on VersaBoard**



The Prone attachment



The Lateral attachment



RT-7045

