

# Multi-Modality Markers



## Catalog Number: MM3007

- Multi-Modality MRI/CT Marker

## Specifications

- 15mm outer diameter
- 3.5mm thick

## Quantity

- 50 markers per jar



## Catalog Number: MM3005

- Multi-Modality Radiology/  
Radiation Therapy Marker

## Specifications

- 15mm outer diameter
- 3.5mm thick
- 5mm central axis hole
- 2mm inner center hole suitable for passage  
of a conventional hypodermic needle

## Quantity

- 50 markers per jar



## Catalog Number: MM3003

- Multi-Modality NucMed/  
PET Marker

## Specifications

- 15mm outer diameter
- 3.5mm thick
- 5mm central axis hole with sealant plastic  
cover that forms a well suitable for the  
injection of a short-life radionuclide  
using a conventional hypodermic needle
- Go to [www.izimed.com](http://www.izimed.com) to see our instructional  
video on how to use the MM3003 marker

## Quantity

- 50 markers per jar



## Multi-Modality Markers

### Benefits of Multi-Modality Markers

- A clear, accurate reference point on any type of scan
- A universal marker
- Shows on CT, MRI, x-Ray, and Mammography
- Shows on all MRI sequences
- Does not indent tissue
- No artifacts
- Disposable

### Benefits Over Other Markers

- Appear on all sequences of MRI even on fat suppression scans
- Does not melt or leak onto the patient
- Size and shape are perfect for extremities and small areas of interest
- No distortions will appear on your images
- Does not indent soft tissue
- Softer design for patient comfort



Quality. Versatility. Innovation.

## Testimonials

“We love using your Multi-Modality Markers because they get in the direct area you want them and do not move around. Yours show up really clear and you know what you are looking at.”

*St. Elizabeth Hospital Health Center*

“The Multi-Modality Marker has better flexibility of use and shows up on all imaging sequences. It stays in the affixed position.”

*Johns Hopkins Hospital*

“Great Markers! Perfect size and shows on both T1 and T2 weighted images. Also, the markers do not indent the skin.”

*Open MRI Specialist of Stockbridge*

“I love the Multi-Modality Markers because they are small and comfortable for the patients. They do not indent soft tissue and they show up great on film.”

*Firelands Community Hospital*

“My radiologist really likes the look of the Multi-Modality Marker and has suggested getting them at other sites. Because of the flat design, they have a much more professional look.”

*Horizon Diagnostic Center*

[www.izimed.com](http://www.izimed.com)